



#### Bibliogrāfija:

- [1] **Parsons M., Newman D.J., Pugia M., Newall R.G., and Price C.P.:** *Performance of a reagent strip device for quantitation of the urine albumin : creatinine ratio in a point of care setting.* Clin Nephrol. 1999 Apr;51(4):220-7.
- [2] **Sasaki M., Pugia M.J., Parker D.R., Kuromoto K., Furukawa I., and Konishi I.:** *Measurement of the albumin content of urinary protein using dipsticks.* J Clin Lab Anal. 1999;13(5):246-50.
- [3] **Kaplan N.M.:** *Microalbuminuria: a risk factor for vascular and renal complications of hypertension.* Am J Med. 1992 Apr 27;92(4B):8S-12S.
- [4] **Apostolovic S., Stanojevic D., Djordjevic V., Jankovic Tomasevic R., Salinger Martinovic S., Djordjevic Radojkovic D., Koracevic G., Stojanovic I., and Pavlovic M.:** *Prognostic significance of microalbuminuria in patients with acute myocardial infarction.* Clin. Lab. 57; 229-235; 2011
- [5] **Pugia M.J., Lott J.A., Clark L.W., Parker D.R., Wallace J.F., Willis T.W.:** *Comparison of urine dipsticks with quantitative methods for microalbuminuria.* Eur J Clin Chem Clin Biochem. 1997 Sep;35(9):693-700.
- [6] **Position Statement: Nephropathy in Diabetes.** Diabetes Care, 27, S79-S83, 2004
- [7] **Pugia M.J., Lott J.A., Profitt J.A., and Cast T.K.:** *High-sensitivity dye binding assay for albumin in urine.* J Clin Lab Anal. 1999;13(4):180-7.
- [8] **Pugia M.J., Lott J.A., Wallace J.F., Cast T.K., Bierbaum L.D.:** *Assay of creatinine using the peroxidase activity of copper-creatinine complexes.* Clin Biochem. 2000 Feb;33(1):63-70.
- [9] **Guy M., Borzomato J.K., Newall R.G., Kalra P.A., and Price C.P.:** *Protein and albumin-to-creatinine ratios in random urines accurately predict 24 h protein and albumin loss in patients with kidney disease.* Ann Clin Biochem. 2009;46: 468-76.
- [10] **Burtis C.A., Ashwood E.R. and Bruns D.E. (eds):** *Tietz Textbook of Clinical Chemistry and Molecular Diagnosis (5th edition).* Elsevier, St. Louis, USA, 2012, 2238 pp, 909 illustrations. ISBN: 978-1-4160-6164-9.



#### Kataloga numurs:

 ACR-9902-1

#### Ražotājs:

  77 Elektronika Kft.  
1116 Budapešta,  
Fehérvári út 98.  
Tel.: (06-1) 206-1480  
Fakss: (06-1) 206-1481  
Epasts: sales@e77.hu  
www.e77.hu

#### Simboli:

-  In vitro diagnostikas medicīnas ierīce
-  Kataloga numurs
-  Partijas numurs
-  Produkts atbilst spēkā esošajām Eiropas Savienības direktīvām.
-  Izlietot līdz datumam
-  Uzglabāšanas temperatūra
-  Ražotājs
-  Sargāt no tiešiem saules stariem vai karstuma
-  Lūdzu, izlasiet iepakojuma instrukciju.
-  Uzmanību
-  Bioloģiskā bīstamība
-  Iepakojumā 25 gabali
-  Tikai vienreizējai lietošanai
-  Nelietot, ja iepakojums ir bojāts
-  Instrukcijas sniegtas angļu valodā
-  Vācu valoda
-  Spāņu valoda